

STANDEZY BED POLE

STRAIGHT HANDLE BED POLE



- Overall Depth 45mm
- Overall Length on frame (upright) 100mm
- Overall Length on frame (folded) 520mm
- Overall Height 690mm
- Handle Height 540mm
- Weight 2.1 kg



The StandEzy Straight Handle Bed Pole allows residents with reduced mobility, but adequate strength, to manoeuvre themselves safely into a sitting position, to the edge of the bed, stand up and transfer into a chair. The stainless steel vertical pole with rubber handle grip allows patients to have confidence that it will provide a strong and stable support.

HOOK HANDLE BED POLE



- Overall Depth 45mm
- Overall Length on frame (upright) 300mm
- Overall Length on frame (folded) 450mm
- Overall Height 580mm
- Handle Height 430mm
- Weight 2.9 kg



The StandEzy Hook Handle Bed Pole provides help to patients with complex mobility issues to get into a sitting position, manoeuvre to the edge of the bed, stand up and transfer into a chair. The stainless steel horizontal handle rotates 360 degrees and provides secure support as the person stands, steps or pivot turns and helps them move with physical support, assisted balance and confidence.

- Maximum Assistance
- Handle has 4 positions
- Clamp or bolt to hospital bed frame
- Safe Working Load 300kg
- Use either side of bed
- Folds away for safety
- Sits inside safety rails
- Raise and lower inside safety rails
- Allows backrest to raise and lower
- Superior strength confirmed by engineering tests
- Chrome plated steel vertical upright pole
- Durable rubber grip

**'STRONG,
STABLE AND
EASY TO OPERATE'**